
















Company	CareFusion		Circadiance		DeVilbiss Healthcare
Product Name	 The Advantage Series Nasal Mask	 SNAPP 2.0	 SleepWeaver ADVANCE Soft Cloth Nasal CPAP Mask with Zzzephyr Seal	 SleepWeaver ADVANCE Soft Cloth Nasal CPAP Mask with new colors and Zzzephyr Seal	 EasyFit Full Face Mask
Website	www.carefusion.com	www.carefusion.com	www.sleepweaver.com	www.sleepweaver.com	www.devilbisshealthcare.com
Price/Warranty	\$47.25/90 days	\$89.50/90 days	\$170 + 30-day money-back guarantee with online sale	\$170 + 30-day money-back guarantee online sale	30-day warranty
Usage					
CPAP	X	X	X	X	X
APAP	X	X	X	X	X
Bi-level	X	X	X	X	X
Size Options	S, M, L	S, M, L, XL, or the starter kit that includes all sizes	One size fits all	One size fits all	S, M, L, XL
Latex Free	X	X	X	X	X
Intentional Leak (ml/s)	None	None	None	None	20 lpm @ 8 cm H ₂ O, 40 lpm @ 20 cm H ₂ O, 50 lpm @ 25 cm H ₂ O
Additional Features		Flexible with no hard plastic on the patient's face or head.	New Zzzephyr Seal, a soft elastic cloth seal designed specifically to fit more comfortably over the patient's nose, now available in three colors.	New Zzzephyr Seal, a soft elastic cloth seal designed specifically to fit more comfortably over the patient's nose, now available in three colors.	Replaceable cushions and forehead pads available in gel or silicone. Exhalation vent is ultra-quiet and designed to direct air parallel to the patient's face and not at a bed partner.
Fitting Support Solutions		Nasal sizing card added with the starter kit.	Call and order DVD or view instructional video online at www.circadiance.com/help	Call and order DVD or view instructional video online at www.circadiance.com/help	Headgear accommodates an array of head profiles. Adjustable forehead support with wider than average forehead cushions for improved customization.
Pressure Range	3-20 cm H ₂ O	3-20 cm H ₂ O	4-20 cm H ₂ O	4-20 cm H ₂ O	4-30 cm H ₂ O
Cleaning Instructions	Wash in warm water with a mild dishwashing soap. Hang to air dry.	Wash in warm water with mild dishwashing soap. Hang and let air dry.	Clean in the clothes washer or hand wash; hang dries quickly.	Clean in the clothes washer or hand wash; hang dries quickly.	Disassemble mask before cleaning. Immerse all components in warm water. Use a mild detergent and soft cloth to remove residue. Clean cavities with a soft brush. Rinse with clean water and allow to air dry.
Headgear Description	The headgear is made with Breathe-O-Prene, a highly breathable material, with Silvertec antimicrobial added to promote patient safety.	The headgear comes in either an adjustable or nonadjustable version and is made with Breathe-O-Prene, a highly breathable material, with Silvertec antimicrobial added to promote patient safety.	One-piece headgear is made entirely of cloth. No plastic buckles, clips, or rigid parts that can scratch and irritate the skin.	One-piece headgear is made entirely of cloth. No plastic buckles, clips, or rigid parts that can scratch and irritate the skin.	The four-point headgear crosses at the back of the patient's head. In contrast to conventional headgear, HEADstrap ensures that no pressure is applied to the patient's neck. It is suitable for every head size.
Interface Description	A nasal mask that offers a triple seal to minimize leaks and a silicone forehead cushion to increase stability and comfort for the patient.	A nasal pillow with cushioning bellows that expand and compress with movement to provide a comfortable reduced leak seal. The two tube design gives the customer the flexibility to sleep on either side without having to readjust after rolling over.	Unlike hard plastic CPAP masks, the SleepWeaver ADVANCE with Zzzephyr Seal inflates like a balloon, softly conforming to the patient's face and creating an intimate seal. Because the elastic cloth is flexible and follows the unique shape of each patient's face, it's more comfortable, doesn't irritate the skin, and doesn't allow air to escape where hard plastic parts don't fit perfectly. Plus, SleepWeaver ADVANCE operates quietly, so that patients (and spouses) can get a good night's sleep.	Unlike hard plastic CPAP masks, the SleepWeaver ADVANCE with Zzzephyr Seal inflates like a balloon, softly conforming to the patient's face and creating an intimate seal. Because the elastic cloth is flexible and follows the unique shape of each patient's face, it's more comfortable, doesn't irritate the skin, and doesn't allow air to escape where hard plastic parts don't fit perfectly. Plus, SleepWeaver ADVANCE operates quietly, so that patients (and spouses) can get a good night's sleep.	The EasyFit Full Face Mask was designed by utilizing computer analyses of facial characteristics from 50,000 individuals globally. Identifying consistencies between most of these facial profiles attributed to the development of the scallop-shaped cushion. This, combined with the lightweight and durable nasal mask, offers a superior fit, seal, and comfort. It also includes the exclusive ball-and-socket joint at the tubing connector, allowing greater freedom of movement, even for the most active sleepers.

Company		DeVilbiss Healthcare	Fisher & Paykel Healthcare		Hans Rudolph Inc	Invacare Corporation
Product Name		 EasyFit Nasal Mask	 Forma Full Face Mask	 Zest Nasal Mask	 V2 Full Face CPAP Mask	 Twilight Full Face Mask
Website		www.devilbisshealthcare.com	www.fphcare.com	www.fphcare.com	www.rudolphkc.com	www.invacare.com
Price/Warranty		30-day warranty	n/a	n/a	6-Month Warranty and Size Guarantee	90-day warranty
Usage	CPAP	X	X	X	X	X
	APAP	X	X	X	X	X
	Bi-level	X	X	X	X	X
Size Options		S, M, L, XL	S, M, L, XL	petite, standard, plus	5 Sizes. L, M, S, XS, and Petite	M
Latex Free		X	X	X	X	X
Intentional Leak (ml/s)		20 lpm @ 8 cm H ₂ O, 40 lpm @ 25 cm H ₂ O	n/a	n/a	n/a	
Additional Features		Replaceable cushions and forehead pads available in gel or silicone. Exhalation vent is quiet and designed to direct air parallel to the patient's face and not at a bed partner.	FlexiFoam Technology	Easy-Clip Silicone Seal with FlexiFoam	Detachable swivel port, airtight seal, new headgear.	Unlimited position of forehead support for customized fit, uses headgear as forehead cushion.
Fitting Support Solutions		Headgear accommodates an array of head profiles. Adjustable forehead support with wider than average forehead cushions for improved customization.	www.vigor8.com	www.vigor8.com	Call or e-mail HRI for pictured fitting instructions.	Adjustable headgear, 360-degree swivel, adjustable forehead support.
Pressure Range		4-30 cm H ₂ O	3 to 20 cm H ₂ O	3 to 25 cm H ₂ O	3 to 40 cm H ₂ O	3-20 cm H ₂ O
Cleaning Instructions		Disassemble mask before cleaning. Immerse all components in warm water. Use a mild detergent and soft cloth to remove residue. Clean cavities with a soft brush. Rinse with clean water and allow to air dry.	www.vigor8.com	www.vigor8.com	Mask, swivel dishwasher safe. Can hand-clean headgear. V2: steam autoclaved, liquid chemical, hot water pasteurized in a hospital setting.	Clean regularly with soap and water.
Headgear Description		The four-point headgear crosses at the back of the patient's head. In contrast to conventional headgear, HEADstrap ensures that no pressure is applied to the patient's neck. It is suitable for every head size and is easy to use.	Stretchgear Headgear is a comfortable and breathable elastic headgear, specially designed to maintain a firm seal.	Stretchgear Headgear is a comfortable and breathable elastic headgear, specially designed to maintain a firm seal.	Simple Strap Foam Headgear with side adjustments and top of the head Tri-Glide adjustment.	Breathable material with soft hook and loop fasteners for adjustability and comfort; quick clips allow the mask to be removed without changing headgear settings.
Interface Description		The EasyFit Nasal Mask was designed by utilizing computer analyses of facial characteristics from 50,000 individuals globally. Identifying consistencies between most of these facial profiles attributed to the development of the scallop-shaped cushion. The mask includes a ball-and-socket joint at the tubing connector, allowing greater freedom of movement.	The Forma Full Face Mask forms a superior seal and provides even more comfort with a greater range of movement through the active contouring of its new FlexiFoam Cushion. A new T-Piece adds stability to the stylized mask base, and works in harmony with the silicone seal and under-chin design to enhance sleep performance.	Small, quiet, and effortless to fit and use, the Zest Nasal Mask offers superb comfort. The innovative design of the Zest Nasal Mask combines the Easy-Clip Silicone Seal with the new FlexiFoam and Advanced Air Diffuser with proven technology.	One piece of molded silicone. No painful forehead support. No hard plastic frame. Goes below the chin to the top of the nose.	Features a simple design and is approved for multipatient use. Dual-walled silicone cushion is easy to remove and clean. By distributing instead of concentrating pressure, the forehead support minimizes telltale red marks and patient complaints. Mask includes frame with adjustable forehead support, cushion, 4 quick clips, and breathable headgear.

Company		Philips Respironics		ResMed		SleepNet Corp
Product Name		 ComfortGel	 FullLife	 Mirage SoftGel	 Swift LT for Her	 MiniMe Pediatric Mask
Website		http://comfortgel.respironics.com	http://fulllife.respironics.com	www.resmed.com	www.resmed.com	www.sleep-net.com
Price/Warranty		90-day warranty for defects of workmanship and materials	90-day warranty for defects of workmanship and materials	90-day warranty	90-day warranty	3-month warranty for manufacturing defects
Usage	CPAP	X	X	X	X	X
	APAP	X	X	X	X	X
	Bi-level	X	X	X	X	X
Size Options		P, S, M, L	S, M, L	M, L	3 pillow sizes: XS, S, M	Single size
Latex Free		X	X	X	X	X
Intentional Leak (ml/s)		approximately 14-50, dependent on PAP pressure (4-30 cm H ₂ O)	approximately 14-50, dependent on PAP pressure (4-30 cm H ₂ O)			
Additional Features		FitPack and DuoPack options available	FitPack and DuoPack options available	Part of ResMed's ConvertAble Series, offering the choice of two mask cushions (Mirage SoftGel and Mirage Activa LT) on one durable frame.	First mask specially designed for women; lightweight (2.3 oz); whisper-quiet comfort with quiet vent design that disperses air gently away from the patient and bed partner.	Mask shell bends to follow facial contours; semirigid tube supports mask position.
Fitting Support Solutions		Fitting guide, video, sizing gauge	Fitting guide, sizing gauge, Web site video clips	Fitting tips online, mask fitting template, and user guide.	Fitting tips online, mask fitting template, and user guide.	Fitting tips and trouble-shooting guide on Web site.
Pressure Range		4-30 cm H ₂ O	4-30 cm H ₂ O	4-20 cm H ₂ O	4-20 cm H ₂ O	3-19 cm H ₂ O
Cleaning Instructions		Mild soapy water without perfumes or lotions. Multiple-patient use lab reprocessing instructions also available.	Mild soapy water without perfumes or lotions. Multiple-patient use lab reprocessing instructions also available.	Hand wash with warm, soapy water (do not use soaps with added scents), rinse thoroughly, and air dry out of direct sunlight.	Hand wash with warm, soapy water (do not use soaps with added scents), rinse thoroughly, and air dry out of direct sunlight.	Hand wash with mild soapy water, rinse thoroughly, and air dry. Multiple-patient use instructions available.
Headgear Description		Headgear with EZ Peel tabs adjusts in four places for a secure, custom fit.	Step-by-step numbered headgear system for easy setup and use; chin support band ensures a secure fit and mask stability; soft headgear has no buckles and keeps the forehead and vision areas open.	ResMed standard Breathe-O-Prene headgear, 3 sizes available.	Breathe-O-Prene headgear with adjustable backstrap can be worn over or under the hair; optional tube retainer allows for side sleeping; expandable stability arms with soft wraps accommodate various facial widths to minimize facial marks and increase comfort.	Cap headgear in S, M, and L.
Interface Description		Nasal mask with replaceable and customizable gel cushion includes adjustable forehead arm, replaceable silicone flap, ball-and-socket headgear attachments. Made of high-strength polycarbonate material with new robust socket wall design.	Minimal contact, full-face mask that is compact and offers a clear line of sight. Dual-textured, grooved cushion creates a robust seal—even with movement. One faceplate fits all sizes.	The DoubleGel cushion's inner blue layer maintains stability and support, while the clear outer layer automatically conforms to the nose and unique facial features of each individual, creating a soft and gentle fit.	The dual-wall nasal pillows seal softly for improved stability and comfort, the rotating barrel allows patient to customize the most comfortable seal, and the small mask frame makes it ideal for side sleeping.	Flexible, bendable mask shell with SoftGel technology.

Information for this guide based on data submitted by product manufacturers. Sleep Review strives for accuracy in all data but cannot be held responsible for claims made by manufacturers.