





















Company		AG Industries		CareFusion	Circadiance	
Product Name		 AG Industries Full Face Mask	 Sopora Nasal Mask	 SNAPP 2.0 Direct Nasal Mask	 SleepWeaver ADVANCE Soft Cloth Nasal CPAP Mask	 SleepWeaver élan Soft Cloth Nasal CPAP Mask
Website		www.agindustries.com	www.agindustries.com	www.carefusion.com	www.circadiance.com	www.circadiance.com
Price/Warranty		\$49.95/30 days	\$39.95/30 days	Masks and headgear free from defects in materials and workmanship for 90 days.	\$170 retail/30-day money-back guarantee	\$170 retail/30-day moneyback guarantee
Usage	CPAP	Yes	Yes	Yes	Yes	Yes
	APAP	Yes	Yes	Yes	Yes	Yes
	Bi-level	Yes	Yes	Yes	Yes	Yes
Size Options		S, M, L	S, M, L	S, M, L, XL	Single size	Regular, Large
Latex Free		Yes	Yes	Yes	Yes	Yes
Intentional Leak		488 ml/s (or 29.3 l/min) @ 10 cm H ₂ O	468 ml/s (or 28.1 l/min) @ 10 cm H ₂ O	30 l/min @ 10 cm H ₂ O	25 l/min @ 10 cm H ₂ O	33 l/min @ 10 cm H ₂ O
Additional Features		Ergonomically designed, odorless, soft, and comfortable. This mask uses CIM technology to create a flexible, self-adjusting forehead bridge for superior comfort and seal.	Ergonomically designed, odorless, soft, and comfortable. This mask uses CIM technology to create a flexible, self-adjusting forehead bridge for superior comfort and seal.	Bellowed nasal interface cushions for patient comfort. Mask cushion flanges provide mask stability and maintain nasal seal for restless sleepers. Conveniently packaged in kit that includes 4 sized cushions.	SleepWeaver lets you fit your patients quickly and comfortably to attain a reliable reading during titration.	All the benefits of the SleepWeaver Advance plus ability to route the tube over the head, nothing in the field of view, replaceable cushion.
Fitting Support Solutions		Fitting instructions are included on mask insert.	Fitting instructions are included on mask insert.	Fitting gauge included in package along with illustrations for fitting tips.	Call and order DVD or view instructional video at www.sleepweaver.com .	All cushion sizes included in starter kit.
Pressure Range		4~20 cm H ₂ O	4~20 cm H ₂ O	4-20 cm H ₂ O	1-20 cm H ₂ O	4-20 cm H ₂ O
Cleaning Instructions		Daily: wipe cushions with wet cloth. Weekly: disassemble components, immerse in warm water/soap solution. Rinse and dry before use.	Daily: wipe cushions with wet cloth. Weekly: disassemble components, immerse in warm water/soap solution. Rinse and dry before use.	Wipe cushions daily with clean wet cloth. Wash weekly, immersing mask and headgear into mild dishwashing soap solution. Dry before use.	Use cold or warm water. Can be washed in clothes washer with detergent or by hand with laundry detergent (always use oil-free products). Hang to dry.	Clean by hand washing with mild nonabrasive soap in warm water. Rinse well. Hang to dry.
Headgear Description		Headgear is soft and comfortable, and features "quick snap" design that makes wearing and cleaning the interface easy.	Headgear is soft and comfortable, and features "quick snap" design that makes wearing and cleaning the interface easy.	Two-point adjustable headgear. Made from AccuMed Breath-O-Prene fabric with SilverTec lining to prohibit bacterial growth.	Headgear is made of flexible Velcro and elastic straps that allow the SleepWeaver mask to be quickly sized to any face.	Four-point adjustment headgear; color coded tabs for top and bottom of headgear.
Interface Description		This mask offers a comfortable, ergonomic design for superior seal. Among the many features that set this mask apart are the flexible bridge, providing comfort and improving patient compliance. Cushion and forehead pads are made of comfortable and odorless medical grade silicone, and headgear is designed for easy wash-and-wear.	This mask offers a comfortable, ergonomic design for superior seal. Among the many features that set this mask apart are the flexible bridge, providing comfort and improving patient compliance. Cushion and forehead pads are made of comfortable and odorless medical grade silicone, and headgear is designed for easy wash-and-wear.	SNAPP 2.0 can be worn with tubing down or up overhead. Nasal cushions form comfortable seal at nares. Quiet exhalation port directs airflow down and away from bed partner. Headgear has 3 adjustment selections, customizing the fit to individual user preference.	The SleepWeaver achieves a good seal and quiet operation using the properties of a balloon. When a balloon is inflated, all points inside the balloon are at the same pressure. This allows the surface areas of the mask to conform to the face, like a balloon, creating the seal.	The SleepWeaver achieves a good seal and quiet operation using the properties of a balloon. When a balloon is inflated, all points inside the balloon are at the same pressure. This allows the surface areas of the mask to conform to the face, like a balloon, creating the seal.

DeVilbiss Healthcare		Fisher & Paykel Healthcare		Global Medical Holdings LLC
				
EasyFit SilkGel Full Face Mask	EasyFit SilkGel Nasal Mask	Pilauro	ZEST Q	ComfortFit Nasal Deluxe+ Mask (100ND+ Series)
www.DeVilbissHealthcare.com	www.DeVilbissHealthcare.com	www.fphcare.com	www.fphcare.com	www.gmhiowa.com
90-day warranty	90-day warranty	N/A	N/A	N/A
Yes	Yes	Yes	Yes	Yes
Yes	Yes	Yes	Yes	Yes
Yes	Yes	Yes	Yes	Yes
S, M, L	S, M, L	One size	Petite, Standard, Plus	S, M, L
Yes	Yes	Yes	Yes	Yes
22.4 lpm @ 8 cm H ₂ O, 37 lpm @ 18 cm H ₂ O	22.4 lpm @ 8 cm H ₂ O, 37 lpm @ 18 cm H ₂ O	28 l/min at 8 cm of pressure	31 l/min at 9 cm of pressure	N/A
Features a soft coating, which provides silky smooth mask and forehead cushions for enhanced comfort. Replaceable cushions and forehead pads available. Exhalation vent is quiet and designed to direct air parallel to the patient's face and not at a bed partner.	Features an ultra-soft coating, which provides silky smooth mask and forehead cushions for enhanced comfort. Replaceable cushions and forehead pads available. Exhalation vent is quiet and designed to direct air parallel to the patient's face and not at a bed partner.	Lightweight at 1.83 ounces; includes flexible tube designed to minimize pressure or pull on the face, and a ball and socket elbow to provide greater flexibility for tube positioning.	The refined Air Diffuser design ensures minimal noise and efficient exhaust dispersal. Use with the cover for extra quiet or without, for a maintenance-free option.	One-piece molded silicone rubber; small, light, and comfortable. Clear mask and frame allow for a good view of the nasal area.
Redesigned headgear accommodates an array of head profiles. Adjustable forehead support with wider than average forehead cushions for improved customization.	Redesigned headgear accommodates an array of head profiles. Adjustable forehead support with wider than average forehead cushions for improved customization.	One convenient size. No manual adjustments necessary. Easy fitting.	Mask fitting guide is provided with every product and is also available for download at www.fphcare.com .	Includes step-by-step instructions that help the patient assemble the mask for maximum fit and comfort.
4-40 cm H ₂ O	4-40 cm H ₂ O	4-20 cm H ₂ O	4-20 cm H ₂ O	4-20 cm H ₂ O
Disassemble mask before cleaning. Immerse all components in warm water. Use a mild detergent and soft cloth to remove residue. Clean cavities with a soft brush. Rinse with clean water and allow to air dry.	Disassemble mask before cleaning. Immerse all components in warm water. Use a mild detergent and soft cloth to remove residue. Clean cavities with a soft brush. Rinse with clean water and allow to air dry.	Wash mask base assembly and silicone seal in soap dissolved in lukewarm water. Rinse thoroughly with fresh water. Leave all parts to dry out of direct sunlight before assembling.	Rinse thoroughly with fresh water; ensure all soap is removed. Dry in direct sunlight before reassembling.	Gently hand wash mask with mild soap and warm water. Rinse mask and components thoroughly. Air dry before reusing.
The four-point headgear crosses at the back of the patient's head. In contrast to conventional headgear, HEADstrap ensures that no pressure is applied to the patient's neck. It is suitable for every head size and is easy to use.	The four-point headgear crosses at the back of the patient's head. In contrast to conventional headgear, HEADstrap ensures that no pressure is applied to the patient's neck. It is suitable for every head size and is easy to use.	The Stretchwise Headgear is made from a soft latex-free elasticized thread with intelligent-stretch properties, and the headgear automatically customizes to the wearer. No manual adjustments necessary.	Stretchgear Headgear Deluxe mask suspension system. A comfortable and breathable elastic headgear, specially designed to maintain a firm seal.	Adjustable, soft headgear with fasteners. Extra buckles provide versatility and movement.
The EasyFit SilkGel Full Face Mask was designed by utilizing computer analyses of facial characteristics from 50,000 individuals globally. Identifying consistencies between most of these facial profiles contributed to the development of the scallop-shaped cushion. The mask includes a ball-and-socket joint at the tubing connector, allowing greater freedom of movement.	The EasyFit Nasal Mask was designed by utilizing computer analyses of facial characteristics from 50,000 individuals globally. Identifying consistencies between most of these facial profiles contributed to the development of the scallop-shaped cushion. The mask includes a ball-and-socket joint at the tubing connector, allowing greater freedom of movement.	New AirPillow Seal Technology: the AirPillow Seal is 3D contoured with fine Variable Thickness Silicone. It uses CPAP pressure to self-inflate to create an effective yet gentle seal into and around the nose. The AirPillow Seal can move and flex in any direction without compromising the seal.	The Zest Q minimizes partner disturbance and ensures a better SleepLife. The F&P Zest Q includes Fisher & Paykel Healthcare features and benefits, and, with the Advanced Air Diffuser, is quieter than ever.	Lightweight, soft, and comfortable. The 100ND+ Series offers forehead support and replacement cushions, maximizing the life of the mask.

Company		Hans Rudolph Inc		InnoMed Technologies	Philips Respironics	
Product Name		 V2 Full Face CPAP Mask	 7500 Series NIV V2 Mask	 Nasal-Aire II	 GoLife for Men GoLife for Women	 TrueBlue
Website		www.rudolphkc.com	www.rudolphkc.com	www.innomedinc.com	www.philips.com/GoLife	www.philips.com/TrueBlue
Price/Warranty		6-month warranty and size guarantee	\$215; 12-month warranty	90-day warranty for defects of workmanship and materials	90-day warranty for defects of workmanship and materials	30-Day Mask Satisfaction Program
Usage	CPAP	Yes	Yes	Yes	Yes	Yes
	APAP	Yes	Yes	Yes	Yes	Yes
	Bi-level	Yes	Yes	Yes	Yes	Yes
Size Options		P, XS, S, M, L	P, XS, S, M, L	XS, S, M, M+, L, XL	P, S, M, L	P, S, M, MW, L
Latex Free		Yes	Yes	Yes	Yes	Yes
Intentional Leak		N/A	N/A	27 l/min @ 10 cm H ₂ O average	4-30 standard liters per minute depending on PAP pressure.	4-30 standard liters per minute depending on PAP pressure
Additional Features		Detachable swivel port, airtight seal, new headgear.	Fits under chin, single piece soft silicone rubber, foam strap type headgear, 5 sizes, leak free.	Features patented L-Flow technology to deliver flow to patient with less turbulence. Ultra-soft nasal inserts.	GoLife cushions are interchangeable. FitPack has all cushions in bag, each wrapped separately.	Freeform spring gives patients the ability to move in any direction. The soft silicone spring helps patients to get a better night's sleep by limiting disturbances to the seal.
Fitting Support Solutions		Call or e-mail HRI for pictured fitting instructions.	Mask sizing gauge available. Call or e-mail HRI for pictured fitting instructions.	Available in an All Sizes Kit. Sizing gauges available. Online videos available. For more information, call (866) 216-7337.	Stand-alone sizing gauge, Web site video clips.	Auto Seal technology is designed to deliver a higher degree of comfort, stability, and freedom of movement with minimal adjustments.
Pressure Range		3-40 cm H ₂ O	0-40 cm H ₂ O	4-20 cm H ₂ O	4-30 cm H ₂ O	N/A
Cleaning Instructions		Mask, swivel dishwasher safe. Can hand-clean headgear. V2: steam autoclaved, liquid chemical, hot water pasteurized in a hospital setting.	Minimal simple two piece disassembly. Dishwasher safe; hot or cold water and soap. Steam autoclaved, hot water pasteurized, liquid chemicals.	Wash daily with warm and mild soapy water. Rinse well. Air dry out of direct sunlight.	Mild soapy water without perfumes or lotions. Multiple-patient use lab reprocessing instructions also available.	Cleaning the mask every day is an important part of maintenance. Wash the components in warm water with mild soap and air dry.
Headgear Description		Simple Strap Foam Headgear with side adjustments and top of the head Tri-Glide adjustment.	3 sizes, soft, lightweight, simple strap type foam headgear.	Features fully adjustable headgear that allows product to be worn in multiple ways, including over the head.	Simple, preformed headgear. Peel tabs easily adjust in four places for a secure, custom fit.	Headgear has inlaid adjustable straps that converge into a single back panel with a crown strap for added stability.
Interface Description		One piece of molded silicone. No painful forehead support. No hard plastic frame. Goes below the chin to the top of the nose.	Anatomically shaped face pieces in 5 different sizes that fit into face and under chin. No hard frame, no forehead cushion, patient can wear glasses and read and watch television.	Cannula style interface with patented sealing within the nose. Virtually eliminates pressure points on the face and bottom of the nose. Comes in an All Sizes Kit or in individual sizes, including a petite sized product line. Replacement frame systems available.	GoLife's soft facial contour arms gently hug the face to maintain stability and a secure seal even as the patient moves. Self-adjusting, angled nasal pillows allow a one-step fit.	TrueBlue—the first blue gel nasal mask with Auto Seal technology. Premium blue gel works together with the thin, form-fitting outer silicone membrane to create an effective, self-adjusting seal.

Company		ResMed		RespCare		Sleepnet
Product Name		 Mirage FX nasal mask	 Quattro FX full face mask	 Aloha Nasal Pillows System	 Hybrid Dual-Airway Mask	 iQ
Website		www.resmed.com	www.resmed.com	www.respcare.com	www.respcare.com	www.sleepnetmasks.com
Price/Warranty		90-day warranty	90-day warranty	90-day warranty for defects of workmanship and materials	90-day warranty for defects of workmanship and materials	90-day limited warranty for manufacturer's defects and a 60-day comfort guarantee.
Usage	CPAP	Yes	Yes	Yes	Yes	Yes
	APAP	Yes	Yes	Yes	Yes	Yes
	Bi-level	Yes	Yes	Yes	Yes	Yes
Size Options		Standard size fits most. Also available in Mirage FX Wide and Mirage FX for Her.	S, M, L	S, M, L	S, M, L	One size
Latex Free		Yes	Yes	Yes	Yes	Yes
Intentional Leak		31 l/min @ 10 cm of pressure	37 l/min @ 10 cm of pressure	28 l/min @ 10 cm H ₂ O	37 l/min @ 10 cm H ₂ O	12-32 standard liters per minute depending on PAP pressure.
Additional Features		Squeeze-tab elbow makes it easy to attach and detach the elbow from the frame. Flexible, "butterfly wing" forehead support allows a dynamic fit range without contacting the skin.	Spring Air cushion and spring frame eliminate need for forehead support, creating an open field of vision. Fewer points of adjustment make fitting fast and straightforward.	Exclusive Arced-Track technology allows the angle of the pillows reservoir to be adjusted to a more customized, natural fit.	Oral cushion portion has a built-in chin flap and a double membrane seal. Pillows portion can be placed in a tall or short position.	Custom Fit Technology for a hand moldable mask. Patented Airgel is soft to the touch and rests directly on the skin for optimal comfort.
Fitting Support Solutions		Mask fitting guides available for download at www.resmed.com. Clinical education services, including mask-fitting workshops. For more information, call (800) 424-0737.	Mask fitting guides available for download at www.resmed.com. Clinical education services, including mask-fitting workshops. For more information, call (800) 424-0737.	Comes in an All Sizes Kit. Includes step-by-step fitting instructions. For more information, call (866) 216-7337.	Comes in an All Sizes Kit. Includes step-by-step fitting instructions and video. For more information, call (866) 216-7337.	Step-by-step instructions with diagrams and instructional videos available on the company website or by calling (800) 742-3646.
Pressure Range		4-20 cm H ₂ O	4-30 cm H ₂ O	4-20 cm H ₂ O	4-30 cm H ₂ O	4-20 cm H ₂ O
Cleaning Instructions		Wash mask system daily with warm, soapy water. Rinse well. Air dry out of direct sunlight. Do not use scented soaps.	Wash mask system daily with warm, soapy water. Rinse well. Air dry out of direct sunlight. Do not use scented soaps.	Wash daily with warm and mild soapy water. Rinse well. Air dry out of direct sunlight.	Wash daily with warm and mild soapy water. Rinse well. Air dry out of direct sunlight.	Wash mask and components in warm water using mild soap. Rinse mask thoroughly. Disinfection method available for multipatient use.
Headgear Description		New, 3D-contoured formfitting SoftEdge headgear is kind to the face and springs to life, ready to fit and wear.	New contoured headgear combined with the spring frame provides stability, durability, and coverage for comfortable, uninterrupted sleep.	Comes standard with Comfort Wraps. Features easy push on/pull out ends for connecting to pillow reservoir.	Five-point adjustable headgear provides for a better seal and increased support over four-point headgears. Features overhead strap for controlled adjustment of pillow seal.	The 3-Point Headgear is a minimalistic design perfect for those using lower pressures. The StableFit Headgear features five adjustment points for a perfect fit.
Interface Description		The Mirage FX is effortless and effective, for comfort through the night. Includes compact frame and soft headgear. With four parts including headgear, wear and care have never been easier.	Full face masks don't have to block a patient's view. Give your patients a mask with a clear field of vision. With a lightweight frame, the Quattro FX offers the coverage and stability of a full face mask in an unobtrusive design.	Full face masks don't have to block a patient's view. Give your patients a mask with a clear field of vision. With a lightweight frame, the Quattro FX offers the coverage and stability of a full face mask in an unobtrusive design.	The Hybrid delivers air to both the nose and the mouth without the typical discomforts of full face masks, such as nasal bridge sores, air leaks near the eyes, and claustrophobic feelings. The ultra-compact design offers clear line of sight with no forehead contact.	Find the perfect fit with Sleepnet's Custom Fit Technology, creating a semi-customizable nasal mask. Patented Airgel is soft to the touch for comfortable therapy. Choose from two headgear designs for more fitting options. The iQ is a high-quality nasal mask made in the USA.