

CATEGORY SPONSORED BY PHILIPS RESPIRONICS

For more than three decades, Philips Respiroics has been a leading global provider of innovative solutions for the sleep and respiratory markets. Its products help manage a broad spectrum of sleep-disordered breathing, chronic obstructive pulmonary disease, asthma, and restrictive lung disorders. From awareness and diagnosis to long-term success, our innovations are designed to enable the best possible outcomes. Solutions that effectively address today's rapidly changing industry landscape and lead to caregiver success require more than just advanced technology. That's why Philips Respiroics has made a commitment to align with caregivers and provide integrated programs and solutions that lead to healthier patients and healthier businesses.

AIRWAY INTERFACES

Company	3B Medical Inc	Airway Management Inc	CareFusion		Circadiance
Product	 iVolve MiniNasal	 myTAP PAP Nasal Pillow Mask <i>See our ad on page 7</i>	 Advantage Series Full Face Mask <i>See our ad on page 3</i>	 SNAPP 2.0 Direct Nasal Mask	 SleepWeaver ADVANCE Soft Cloth Nasal CPAP Mask
Website	www.3Bproducts.com	www.mytappap.com	www.carefusion.com	www.carefusion.com	www.circadiance.com
Warranty (months)	3	3 (patient compliance guarantee)	1	3	3
Available Sizes	S, M, L	S, M, L nasal pillows included	S, M, L, XL	S, M, L, XL	S, one size fits most
Latex-free	X	X	X	X	X
Intentional Leak @ 10 cm H₂O (l/min)	29	30	Info not available	30	25
Fitting Support Solutions	Comes with instructional booklet with fitting instructions and troubleshooting guide.	Patient and clinician fitting guide included; fitting video on myTAP-PAP.com and DVD available.	Instructions for use included with purchase.	Fitting gauge included in package along with illustrations for fitting tips.	www.circadiance.com ; DVD is included with the mask.
Pressure Range (cm H₂O)	4-30 (recommended)	4-18	3-20	4-20	4-20
Cleaning Instructions	Hand wash by gently scrubbing in warm water using mild soap. Use a soft-bristled brush to clean the port cover. Rinse all components well with water and allow to air dry out of direct sunlight.	After each use, rinse the mouth-piece in cold water and wipe the pillows' outer portion with a damp cloth. Once a week, disassemble the mask and soak parts for no longer than 10 minutes in cold water and mild liquid soap. Remove and rinse in cold water then allow the parts to air dry.	Clean with mild detergent. It does not require disassembly.	Wipe cushions daily with clean wet cloth. Wash weekly, immersing mask and headgear into mild dishwashing soap solution. Dry before use.	Clean by hand washing with mild nonabrasive soap in warm water. Rinse well. Hang to dry.
Headgear Description	Quick-release headgear with clips that simplify attachment and release.	Headgear is available for the myTAP PAP Nasal Pillow Mask. However, the mask does stay securely in place without straps.	Adjustable headgear.	Two-point adjustable headgear. Made from AccuMed Breath O-Prene fabric with SilverTec lining to prohibit bacterial growth.	Headgear is soft, comfortable, and features a 3-point adjustable strap system.
Interface Description	Composed of only 2 parts, the mask is extremely easy to fit. Its lightweight and small profile makes it a perfect mask for patients looking for a lighter and less obtrusive mask.	The myTAP PAP Nasal Pillow Mask is a strapless nasal pillow mask that is intended for individuals who have been diagnosed as requiring CPAP or bi-level ventilator treatment. The myTAP PAP is intended for single patient adult use in the home, hospital, or other clinical setting.	Seal around the nose and mouth with triple layer seal design.	SNAPP 2.0 can be worn with tubing down or up overhead. Nasal cushions form comfortable seal at nares. Quiet exhalation port directs airflow down and away from bed partner. Headgear has 3 adjustment selections, customizing the fit to individual user preference.	The SleepWeaver uses the properties of a balloon. When a balloon is inflated, all points inside the balloon are at the same pressure. This allows the surface areas of the mask to conform to the face.
Additional Features	The iVolve MiniNasal has quiet venting and very soft silicone cushioning to improve patient comfort.	The myTAP PAP Nasal Pillow Mask encourages patient compliance by solving the problems commonly associated with traditional CPAP Masks—from awkward straps to disruptive air leaks. The unique and innovative design of the myTAP PAP NPM delivers unparalleled performance and is backed by a Compliance Guarantee Program.	Elbow swivel 360 degrees with anti-asphyxia valve, built-in vent openings.	Bellowed nasal interface cushions for patient comfort. Mask cushion flanges provide mask stability and maintain nasal seal for restless sleepers. Conveniently packaged in kit that includes 4 sized cushions.	SleepWeaver lets you fit your patients quickly and comfortably to attain a reliable reading during titration.

Information for this guide based on data submitted by product manufacturers. Sleep Review strives for accuracy in all data but cannot be held responsible for claims made by manufacturers.

MATRIX SPONSOR

Fisher & Paykel Healthcare	Hans Rudolph Inc		InnoMed Technologies	Philips Respironics	ResMed	
	 See our ad on page 15				 See our ad on page 23	
Simplus Full Face Mask	V2 Full Face CPAP Mask	QUEST HOMECARE Full Face Mask	Aloha	Nuance and Nuance Pro, gel pillows mask	AirFit P10	AirFit F10
www.fphcare.com/simplus	www.rudolphkc.com	www.rudolphkc.com	www.innomedinc.com	www.theanswer.philips.com/masks/nuance	www.resmed.com/us/en/consumer/products/masks/airfit-p10.html	www.resmed.com/us/en/consumer/products/masks/airfit-f10.html
3	6 (includes size guarantee)	6	3	3 (limited warranty)	3	3
S, M, L	P, XS, S, M, L	P, XS, S, M, L	S, M, L pillows included	S, M, L	XS, S, M, L	XS, S, M, L
X	X	X	X	X	X	X
35	N/A	No, pressure tested	~29	21.9	33	37
Fitting gauge included; videos at www.youtube.com/fphcare .	Call or e-mail HRI for pictured fitting instructions.	Mask sizing gauge, fitting instructions included; videos at www.rudolphkc.com .	Adjustable parts. Instructions in manual.	www.theanswer.philips.com/nuance	www.resmed.com/us/en/consumer/support/masks/how-to-fit-your-mask.html	www.resmed.com/us/en/consumer/support/masks/how-to-fit-your-mask.html
4-25	3-40	3-40	3-20	4-20	4-20	3-40
Minimal parts to clean; soap in lukewarm water; high level disinfection available.	Mask, swivel dishwasher safe. Can hand-clean headgear. V2: steam autoclaved, liquid chemical, hot water pasteurized in a hospital setting.	Soak components in warm water with mild detergent, sponge or soft brush, rinse in warm water, air dry, or use lint-free cloth.	See the manual.	Hand wash in warm water with a mild liquid dish washing detergent. Rinse thoroughly. Air dry completely before use. Make sure the entire mask is dry before use. Lay the fabric parts flat or line dry.	Hand wash each day with mild soap; avoid harsh cleaning products; allow to dry away from direct sunlight.	Gently hand wash with mild soap; avoid harsh cleaning products; allow to dry away from direct sunlight.
ErgoForm breathable headgear incorporates stretch and non-stretch panels. It holds its shape and allows stretch where required for optimal usability and performance. Embedded Velcro tabs with easy-hook clip hooks enable fast fitting and removal. Crown strap for stability.	Simple strap foam headgear with side adjustments and top of the head tri-glide adjustment.	Light gray-on-gray color; soft foam and nylon; mounts with simple 4 clips snapped into place; secure comfortable fit.	Reinforced with over-the-head strap. Tubing can be secured in multiple locations.	All fabric non-slip headgear.	Includes soft and stable QuickFit headgear, which auto-adjusts to provide gentle stability and is designed with a split strap for minimal facial contact.	The seamless SoftEdge headgear is gentle on the skin and features premium rolled-edge fabric for reduced facial markings and easy fit. The SoftEdge headgear requires minimal tension to achieve a comfortable fit.
The F&P Simplus incorporates three key components: the RollFit Seal, ErgoForm Headgear, and Easy Frame, all designed to work in harmony. In combination, these components offer the comfort, seal, and easy use to simply fit and perform.	One piece of molded silicone. No painful forehead support. No hard plastic frame. Goes below the chin to the top of the nose.	Fits into the face securely for leak-free fit, won't ride up face due to chin cup, no hard frame so no skin breakdown, 5 sizes, glasses allowed while wearing.	Nasal pillow system.	Philips Respironics brings the comfort and sealing power of gel to nasal pillows. As the latest in our groundbreaking new line of masks, Nuance gel pillows represent the new choice for better sleep therapy for clinicians, home care providers, and patients.	Nasal pillows adapt naturally to movement during sleep, maintaining a secure fit. Its mesh vent, QuietAir, diffuses air with minimal disturbance to patients or their bed partners, making the AirFit P10 50% quieter than its predecessor.	The AirFit F10's four-part design makes it easy to take apart, clean, and reassemble. It is quiet, lightweight, and provides a clear line of sight. The dual wall Spring Air cushion enhances comfort at the nasal bridge and maintains a high performance seal.
Integrated diffuser that minimizes noise and draft for the bed partner, lightweight frame with easy clip, ball and socket swivel connector.	Detachable swivel port, airtight seal, new headgear.	Low cost home care mask allows all patients to afford to wear the high quality Hans CPAP/BiLevel mask. Distributor pricing available.	Exclusive Arced-Track technology adjusts angles for natural fit. Floating ball-and-socket elbow allows various sleeping positions. Incredibly stable, comfortable seal.	Fabric or gel frame, non-slip headgear, gel pillows provide greater comfort and seal, flexible tubing, conforms to different size nostrils.	The ResMed AirFit P10 is light, quiet, comfortable, and easy to use and assemble. It features dual wall nasal pillows for gentle fit and secure seal, minimalist design, and intelligent enhancements for men and women.	AirFit F10 full face mask is comfortable, reliable, and easy to use. It has a sleek design, and is the lightest full face mask on the market. A high performance seal delivers effective therapy.